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Specifications

Table of Contents

	Page
I . Product summary	1
II . Features	1
III . Specifications	2
IV . Functions	8
V . Standard input signal	1 0
VI . External Dimensions	1 1
VII . Attachments	1 2
VIII . Warranty	1 2
● Cases For Indemnity (Limited Warranty)	

TOSHIBA TELI CORPORATION

I . Product summary

This monitor (**T24MSA001-MD**) is employed a 24-inch wide TFT-LCD display panel and Color LCD monitor for medical use. This LCD monitor features high brightness, wide view angle, high resolution (WUXGA) and displaying 16.77 million colors. The LCD monitor may be displayed a single screen or double screen to display VIDEO picture (NTSC, PAL) and PC picture (VGA~UXGA, WUXGA). To control the monitor, in addition to the OSD control from the key switch on the front of the unit, it may also be used external communication (RS-232C) for control. The power supply used for the monitor is designated AC adaptor that meets medical standards.

II . Features

1 . Accepting NTSC/PAL's TV signal in VIDEO

The built in scaling function allows the monitor to display 4:3 aspect ratio NTSC/PAL pictures into full screen or fixed aspect ratio.

The monitor shows either JUST/UNDER/OVER SCAN at fixed aspect ratio.

2 . Accepting VGA ~ WUXGA signal VIDEO PC

The monitor automatically detects the size of picture of the input signal and displays it correctly.

The built in scaling function allows the monitor to convert VGA, SVGA, XGA, SXGA or UXGA to WUXGA equivalent and display it in WUXGA, or display in fixed aspect ratio or equal ratio magnification.

3 . Dual screen display

In addition to a single screen display, dual screen display, Picture In Picture (PIP) or Picture Out Picture (POP) is available.

4 . Adjustment and control function

Adjustments and control through On Screen Display (OSD) or external control by RS-232C are available.

5 . Color temperature setting

The color temperature selection in 9300K / 6500K / USER is available.

In USR mode, R / G / B independently adjustment is available.

6 . Reduced light brightness change during SLEEP mode reset

When no signal is input to the monitor for preset duration, the monitor goes to sleep mode without turning backlight off. About two hours later, the monitor goes to energy saving mode.

7 . Meet VESA specifications, support free mount

Since screw holes in the back of the unit are VESA standard, in addition to the self-standing stand, a wall mount or arm mount can be used.

※ This monitor supports 100mm x 100mm and 200mm x 100mm pitch.

※There are indication restrictions. Please confirm to instruction manual about detail usage.

III. Specifications

1. LCD display panel portion

- | | | |
|-----------------------------|---|------------------------------------------------------------------------------------------------------------------|
| 1) Display device | : | Model 24 a-Si TFT-LCD module |
| 2) Display method | : | IPS method |
| 3) Display screen size | : | 518.4 (H) × 324.0 (V) [mm] (diagonal equivalent to model 24) |
| 4) Number of displayed dots | : | 1920 (H) × 1200 (V) [Dot] |
| 5) Picture pitch | : | 0.270 (H) × 0.270 (V) [mm] |
| 6) Picture pitch layout | : | RGB vertical stripe arrangement |
| 7) Display color | : | 167000000 color (RGB 8 Bit each) |
| 8) Back light | : | Direct light method 6 cold cathode fluorescent lamp
(Backlight single unit service life 50,000 hours/Ta=25°C) |

Caution : Since the backlight (back light) is dispensable part, the life expectancy is for reference only.

This number is not a guarantee. Since the life of the backlight varies significantly by the temperature where the monitor is placed, please use carefully.

For details, please contact our company's sales department.

- | | | |
|------------------------|---|-------------------------------------------------|
| 9) Contrast | : | 800: 1(Typ) |
| 10) visual field angle | : | Up, Down, Right and Left, 89 degrees each (Typ) |
| 11) Response speed | : | Rise 6.5 ms (Typ)
Fall 7.5 ms (Typ) |

Caution : The numbers listed above are guaranteed by the LCD module manufacture (Ta =25°C).

12) Spot/line defects

		Item	Specifications
1.	Line defect		Nothing
2.	Bright spot		≦ 3 dots
3.	Dark spot		≦ 5 dots
4.	Total spot defects		≦ 8 dots
5.	Continuous spot defect	2 continuous blight spots	1 pair
		2 continuous bark spots	2 pair
		Distance between 2 light spots	≧ 15 mm
		Distance between 2 black spots	≧ 15 mm

Caution : ·Pixel defects (bright spot/dark spot) counts one sub pixel (dot: each RGB) as one pixel.

·Two continuous bright spots are that they form two continuous dots (horizontal, vertical or diagonal).

·(Spot defects) defect zone: when defects exceeding 1/2 of one sub pixel, it is considered as one defect spot to carry out counting.

·Adjacent defects are checked with one of the white, black, green, blue, red and gray patterns to carry out check.

13) Front protect glass

- | | | |
|-------------------|---|---------------|
| Materials | : | Acrylic resin |
| Surface treatment | : | Hard coating |

* The information in this documents is subject to change without notice.

14) Cosmetic appearance specifications

Check Item	Criteria			
Scratch	Width: W (mm)	Number(n)	Length: L (mm)	Number(n)
	$W \leq 0.1$	Disregarded	$L < 0.5$	Disregarded
	$0.1 \leq W \leq 0.2$	$n \leq 9$	$0.5 \leq L < 1.0$	$n \leq 9$
			$1.0 \leq L < 7.0$	$n \leq 4$
Dent Air bubble	Average diameter: D(mm)	Number(n)		
	$D < 0.15$	Disregarded		
	$0.15 \leq D < 0.25$	Disregarded (Distance between defects space it more than 50mm.)		
	$0.25 \leq D \leq 0.5$	$n \leq 9$		
Foreign object (Linear)	Width: W(mm)	Number(n)	Length: L (mm)	Number(n)
	$W < 0.01$	Disregarded	$L < 0.3$	Disregarded
	$0.01 \leq W \leq 0.20$	$n \leq 6$	$0.3 \leq L \leq 4.0$	$n \leq 6$
Foreign object (Circular)	Average diameter: D(mm)	Number(n)		
	$D < 0.15$	disregarded		
	$0.15 \leq D < 0.25$	Disregarded (Distance between defects space it more than 50mm.)		
	$0.25 \leq D \leq 0.5$	$n \leq 9$		
	$0.5 \leq D \leq 0.70$	$n \leq 4$		

Caution : Determination specifications※1 Ambient temperature: $25 \pm 5^{\circ}\text{C}$ ※2 Humidity: $65 \pm 10\%$ RH.

※3 Ambient light brightness: 300~700 Lx

※4 The check must be done without panel protection cover.

※5 When problem occurs the check, re-check the display panel during start up.

Caution : The above specifications are factory outgoing specifications. During use, the bright spots, dark spots, scratch, wrinkles, dent, dust and foreign objects etc. conditions discovered by customers are out of warranty.

2. Image signal input portion

2 - 1. VBS signal (Composite signal)

- | | | | |
|--------------------|--------------|---|---------------------------------------------------------------------------|
| 1) Input number | Image signal | : | VBS 1 system |
| 2) Input Connector | | : | BNC Connector 1 system |
| 3) Input signal | method | : | Meet NTSC signal and PAL signal (Interlace) |
| 4) Input level | Image signal | : | 1. $0 \pm 0.05V_{(p-p)}$
(VIDEO: $0.7V_{(p-p)}$ SYNC: $0.3V_{(p-p)}$) |
- ※ When output connector is not connected, or output connector is terminated with 75Ω .

2 - 2. Y/C signal


- | | | | |
|--------------------|--------------|---|--------------------------------------------------------------------------------------------------------------------|
| 1) Input number | Image signal | : | Y/C 1 system |
| 2) Input Connector | | : | Y/C Connector 1 system |
| 3) Input signal | method | : | Meet NTSC signal and PAL signal (interlace) |
| 4) Input level | Image signal | : | Y 1.0 $\pm 0.05V_{(p-p)}$ (75 Ω internal fixed)
C 0.286 $\pm 0.02V_{(p-p)}$ (75 Ω internal fixed) |

* The information in this documents is subject to change without notice.

2 - 3. Analog RGB signal

- 1) Input number Image signal :
 Synchronize signal :
 2) Input Connector :

RGB 1 system each
 HD·VD 1 system each
 D-sub 15pin (male) 1 system

Front View

Pin	Label	Description
1	R	Red Image signal
2	G	Green Image signal
3	B	Blue Image signal
4	GND	Ground
5	GND	Ground
6	AGND	Analog Ground (R)
7	AGND	Analog Ground (G)
8	AGND	Analog Ground (B)
9	D5V	+5V (for DDC)
10	GND	Ground
11	GND	Ground
12	DDC SDA	DDC serial data year (Serial data)
13	HS	Horizontal Synchronize signal
14	VS	Vertical Synchronize signal
15	DDC SCL	DDC serial type clock (serial clock)


- 3) Input signal Image signal :
 Synchronize signal :
 4) Input (level) Image signal :
 Synchronize signal :
 5) Plug and Play function :

Analogue RGB signal (progressive)
 Vertical drive signal, Horizontal drive signal
 (Separate Sync) positive and negative two polarities
 RGB 0.7 ± 0.05V_(p-p) (75Ω internal fixed)
 HD·VD TTL level (1kΩ internal fixed)
 Meet VESA DDC 2B specifications

2 - 4. DVI signal

- 1) Input number
 Digital signal Image signal :
 Control signal :
 Analogue signal Image signal :
 Synchronize signal :
 2) Input connector :

T.M.D.S. (Single Link) 1 system each
 DDC, Hot Plug Detect 1 system each
 RGB 1 system each
 HD·VD 1 system each
 DVI 29pin (female) 1 system

Front View

Pin	Description	Pin	Description
1	T.M.D.S. Data2-	16	Hot Plug Detect
2	T.M.D.S. Data2+	17	T.M.D.S. Data0-
3	GND	18	T.M.D.S. Data0+
4	OPEN	19	GND
5	OPEN	20	OPEN
6	DDC Clock	21	OPEN
7	DDC Data	22	GND
8	Analog VS	23	T.M.D.S. Clock+
9	T.M.D.S. Data1-	24	T.M.D.S. Clock-
10	T.M.D.S. Data1+	C1	Analog Red
11	GND	C2	Analog Green
12	OPEN	C3	Analog Blue
13	OPEN	C4	Analog HS
14	+5V power	C5	GND
15	GND		

* The information in this documents is subject to change without notice.

3) Input signal

Digital signal	:	Signal transmission method T.M.D.S. (Transition Minimized Differential Signaling)
Analog signal image signal	:	Analog RGB signal (progressive)
Synchronize signal	:	Vertical drive signal, Horizontal drive signal (Separate Sync) positive and negative two polarities

4) Input level:



Digital signal	:	Meet DVI (Digital visual interface specification Rev.1.0)
Analog signal Image signal	:	RGB $0.7 \pm 0.05V_{(p-p)}$ (75 Ω internal fixed)
Synchronize signal	:	HD-VD level (1k Ω internal fixed)

5) Plug and Play function

Meet VESA DDC 2B specifications

2 - 5. Serial communication signal

1) Communication method	:	RS - 232 communication
2) Communication format	:	Baud rate 9600/169200/38400bps(bit per sec.) Data length 8 bit, binary Parity Non / Odd / Even Stop bit 1 bit, binary Code type 8 bit, binary, ASCII type
3) Communication connector	:	D-sub 9pin (female) 1 system

4) Communication cable : straight type

3. Image signal output portion



3 - 1. VBS signal (Composite signal)

1) Output number Image signal	:	VBS 1 system
2) Output connection	:	BNC Connector 1 system
3) Output signal Method	:	Based on image input signal
4) Input level image signal	:	Based on image input signal (internal input, output) ※ When connecting output connectors, open internal 75 Ω terminal: automatic termination

4. Input power portion

- Monitor portion

- | | | |
|---------------------------|---|-------------------------------------------------------------------------------------------------------------------------------|
| 1) Input electric voltage | : | DC + 24V±5% (DC + 22. 8~25. 2V) |
| 2) Current consumption | : | Under 4.6A |
| 3) Ripple noise | : | Under 300 mV _(P-P) |
| 4) Connectors used | : | C4255WR-E-2X06P
(Manufactured by Shanchih Enterprise Company) |
| 5) Connector suitable | : | C4255HF-2X06P
(Manufactured by Shanchih Enterprise Company)
39-01-2060
(Manufactured by MOLEX Company) or equivalent |

- When using accessories AC Adaptor

- | | | |
|---------------------------------------|---|----------------------------------------------------------------------------|
| 1) Input voltage/electric consumption | : | AC90~110V / 120W and under
AC110~264V / 118W and under
※Display only |
| 2) frequency | : | 50/60Hz |
| 3) Phase | : | Single phase |

©Specifications of the AC ADAPTOR

Made by PROTEK Company : PMP150-14-J2

Input	▪ Voltage	:	AC100 ~ 240V ± 10%
	▪ Frequency	:	50 / 60Hz
	▪ Surge current	:	Under 60A (AC115V) Under 120A (AC230V, 25°C Cold start)
	▪ Ground leakage electric current	:	0. 09mA (max) (AC115V、 60Hz on input) 0. 15mA (max) (AC230V、 50Hz on input)
Output	▪ Voltage	:	DC+24V ± 5%
	▪ Minimum load current	:	0. 0A
	▪ Maximum load current	:	6. 25A
	▪ Output range	:	22. 8 ~ 25. 2V
	▪ Output replay-noise	:	480mV _(P-P) (max)
	▪ Power factor	:	Over 85%
	▪ Load regulation	:	± 5% (max)
	▪ Over current protection	:	Theo output is protected to short circuit conditions.
	▪ Over voltage protection	:	When at over voltage protection activates (standard out put force under 112 ~ 140%), it stops outputting.
	▪ DC cable length (Cable)	:	1. 8 ± 0. 05m (DC connector are not included)
	▪ Dimensions	:	207.6±0.5(W)×49.7±0.5 (H)×82.0±0.5 (D)[mm] (rubber foot and extended parts are not included)
	▪ Weight	:	0.97±0.30kg

* The information in this documents is subject to change without notice.

- Caution** : When connecting the monitor to the power, make user to use this AC adaptor.
- Caution** : The DC cable connecting to this product shall use UL certified standard parts that can sustain voltage over DC40V and electric current over 6.0 A, and use temperature over 60°C.
- Caution** : An AC Cable is not included. It is recommend to use cable that meets UL2601-1, CAN/CSA-C22.2NO.601 requirement and IEC60601-1 requirement.
- Check the power line specifications of the are that the monitor is sued.

[Reference]			
Japan	Voltage	Single phase 100 V	Frequency 50/60Hz
USA	Voltage	Single phase 220 V	Frequency 60 Hz
European union(E.U)	Voltage	Single phase 220 V - 240 V	Frequency 50 Hz

5. Operation, storage and shipping conditions

1) Use

- Ambient temperature : 0 ~ +35°C
- Ambient humidity : 20 ~ 85% (without due)
- Ambient pressure : 70kPa ~ 106kPa
- Use condition : Indoor and fixed location use

2) Storage

- Storage temperature : - 10 ~ +60°C
- Storage humidity : 20 ~ 85% (without due)
- Storage pressure : 70kPa ~ 106kPa

3) Shipping

- Shipping temperature : - 10 ~ +60°C
- Shipping humidity : 20 ~ 85% (without due)
- Shipping pressure : 70kPa ~ 106kPa

6. Dimensions

- 1) Monitor : 580.4±3 (W) × 421.3±3 (H) × 83.5±2 (D) mm
(Not include extended parts such as switch and connector etc.)

7. Weight

- 1) Monitor : Around 8. 4 kg
(accessories and cables are not included)

8. Applicable specifications

1) Safety standards

(medical device electric appliance specifications) : IEC60601-1

※Classification :

- Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide
- Continuous operation

2) EMC specifications : IEC60601-1-2

3) Enclosure specifications : IP-X1

※ When install this product perpendicular only, it is applicable to a specification of IP-X1.

4) Certifications : GM, cTUVus, CE

5) FCC (Federal Communications Commission) : FCC sub Part 15 Class B

VESA, DPMS, and DDC are the trademarks of Video Electronics Standards Association.

TMDS is registered trademark of Silicon Image, Inc.

Other company names, and product names are trademarks or registered trademarks of individual companies.

* The information in this documents is subject to change without notice.

IV. Function

1. OSD (On Screen Display) control

To control the monitor, in addition to the OSD control from the key switch on the front of the unit and available all setting While seeing an OSD menu on an LCD screen.

About an input image signal, the following adjustment are available.

- 1) CVBS / S-VIDEO SETUP (Composite signal or Y/C signal input)
 - ① Contrast adjustment ② Brightness adjustment
 - ③ Color adjustment ④ Hue adjustment
 - ⑤ Color temperature setting (9300K / 6500K / USER adjustment)
 - ⑥ Motion compensation level ⑦ Video sharp adjustment
 - ⑧ AGC setting ⑨ Scale mode setting (FULL / UNDER / OVER / JUST)
 - ⑩ Horizontal position of screen display adjustment
 - ⑪ Vertical position of screen display adjustment
 - ⑫ Default (Reset selected SETUP MENU parameters to factory defaults.)
- 2) VGA / DVI-D / DVI-A SETUP (Analog RGB signal or DVI-D signal or DVI-I signal input)
 - ① Auto adjust (Display position auto adjustment function)
 - ② Contrast adjustment ③ Brightness adjustment
 - ④ Color temperature setting (9300K / 6500K / USER adjustment)
 - ⑤ Sharp adjustment ⑥ Phase adjustment
 - ⑦ Horizontal total dot adjustment
 - ⑧ Scale mode setting (FULL / ASPECT / ONE TO ONE)
 - ⑨ Horizontal position of screen display adjustment
 - ⑩ Vertical position of screen display adjustment
 - ⑪ Default (Reset selected SETUP MENU parameters to factory defaults.)
- 3) PIP / POP SETUP (All input signal is common)
 - ① PIP/POP mode setting (PIP / POP / OFF)
 - ② PIP/POP source setting (Set an input signal on PIP/POP screen.)
 - ③ Horizontal position of PIP screen display adjustment
 - ④ Vertical position of PIP screen display adjustment
 - ⑤ SWAP function
 - ⑥ Default (Resets "PIP/POP SETUP" parameters to factory defaults.)

		MAIN				
		CVBS	S-VIDEO	VGA	DVI-Analog	DVI-Digital
SUB	CVBS		○	○	○	○
	S-VIDEO	○		○	○	○
	VGA	○	○		×	×
	DVI-Analog	○	○	×		×
	DVI-Digital	○	○	×	×	

Available signal combinations in PIP mode

		MAIN (LEFT)				
		CVBS	S-VIDEO	VGA	DVI-Analog	DVI-Digital
SUB(RIGHT)	CVBS		○	○	○	○
	S-VIDEO	○		○	○	○
	VGA	○	○		×	×
	DVI-Analog	○	○	×		×
	DVI-Digital	○	○	×	×	

Available signal combination in POP mode

* The information in this documents is subject to change without notice.

4) GENERAL SETUP (All input signal is common)

- ①AUTO SOURCE DETECT setting
- ②Backlight setting
- ③Channel display setting (NORMAL ON / SLIM ON / OFF)
- ④VIDEO COMB FILTER setting
- ⑤Horizontal position of OSD screen display adjustment
- ⑥Vertical position of OSD screen display adjustment
- ⑦OSD display timeout setting
- ⑧External control is available setting
- ⑨External control baud rate setting (9600 / 19200 / 38400)
- ⑩External control parity setting (NON / ODD / EVEN)
- ⑪Default (Resets "GENERAL SETUP" parameters to factory defaults.)
- ⑫All default (Resets all SETUP MENU parameter to factory defaults)

2. Memory function

Memorize a setup item adjusting in OSD or external control.

※Reset item parameters to factory defaults by OSD or external control.

3. External control (RS232C)

Serial communication connectors in the back of the unit, it may also be used external communication by RS-232C (setting by OSD equivalent) are available.

4. Energy saving mode

When no signal is input to the monitor for preset duration, the monitor goes to sleep mode without turning backlight off. About two hours later, the monitor goes to energy saving mode.

5. support free mount


Since screw holes in the back of the unit are VESA FPMPMI standard, in addition to the self-standing stand, a wall mount or arm mount can be used.

※ **This monitor supports 100mm x 100mm and 200mm x 100mm pitch.**


6. Cable prevention of falling off function

The prevent Y/C cable from falling off by using attached Nylon clamp and Bind screw.


【Nylon clamp】

- TOP ANGLE VIEW -

- SIDE ANGLE VIEW -

- DIAGONAL ANGLE VIEW -

※ Please use fixing screws of the individual cables to fix the image use cable and communication cable to this product.
(Please use inch screw (#4-40))

V. Standard input signal

●PC Signal (VGA、DVI-A、DVI-D)

No.	Resolution		Dot CLK (MHz)	Horizontal frequency			Vertical frequency			Polarity		Signal name Type	VGA	DVI-A	DVI-D
	Horizontal	Vertical		H Freq	Sync Time	Front Porch	V Freq	Sync Time	Front Porch	H	V				
	(Dot)	(H)		(KHz)	(Dot)	(Dot)	(Hz)	(H)	(H)						
1	640	400	25.175	31.469	96	16	70.086	2	12	N	N	IBM VGA	○	○	○
2	640	400	31.500	37.861	64	32	85.080	3	1	N	P	VESA Standard	○	○	○
3	640	480	25.175	31.469	96	16	59.940	2	10	N	N	Industry	○	○	○
4	640	480	31.500	37.861	40	24	72.809	3	9	N	N	VESA Standard	○	○	○
5	640	480	31.500	37.500	64	16	75.000	3	1	N	N	VESA Standard	○	○	○
6	768	480	27.500	29.364	32	48	59.987	6	2	P	N	CVT Red.Blanking	○	○	○
7	800	600	36.000	35.156	72	24	56.250	2	1	P	P	VESA Guidelines	○	○	○
8	800	600	38.250	37.354	80	32	59.861	4	3	N	P	CVT	○	○	○
9	800	600	40.000	37.879	128	40	60.317	4	1	P	P	VESA Guidelines	○	○	○
10	800	600	50.000	48.077	120	56	72.188	6	37	P	P	VESA Standard	○	○	○
11	800	600	49.500	46.875	80	16	75.000	3	1	P	P	VESA Standard	○	○	○
12	1024	768	65.000	48.363	136	24	60.004	6	3	N	N	VESA Guidelines	○	○	○
13	1024	768	75.000	56.476	136	24	70.069	6	3	N	N	VESA Standard	○	○	○
14	1152	864	108.000	67.500	128	64	75.000	3	1	P	P	VESA Standard	○	○	○
15	1280	720	74.500	44.772	128	64	59.855	5	3	N	P	CVT	○	○	○
16	1280	720	95.750	56.456	128	80	74.777	5	3	N	P	CVT	○	○	○
17	1280	768	68.250	47.396	32	48	59.995	7	2	P	N	CVT Red.Blanking	○	○	○
18	1280	768	117.500	68.633	136	80	84.837	7	3	N	P	CVT	○	○	○
19	1280	960	101.250	59.699	128	80	59.939	4	3	N	P	CVT	○	○	○
20	1280	960	108.000	60.000	112	96	60.000	3	1	P	P	VESA Standard	○	○	○
21	1280	960	148.500	85.938	160	64	85.002	3	1	P	P	VESA Standard	○	○	○
22	1280	1024	109.000	63.668	136	80	59.895	7	3	N	P	CVT	○	○	○
23	1280	1024	108.000	63.981	112	48	60.020	3	1	P	P	VESA Standard	○	○	○
24	1280	1024	135.000	79.976	144	16	75.025	3	1	P	P	VESA Standard	○	○	○
25	1280	1024	157.500	91.146	160	64	85.024	3	1	P	P	VESA Standard	○	○	○
26	1360	768	84.750	47.720	136	72	59.799	5	3	N	P	CVT	○	○	○
27	1360	768	109.000	60.288	144	80	74.891	5	3	N	P	CVT	○	○	○
28	1400	1050	121.750	65.317	144	88	59.978	4	3	N	P	CVT	○	○	○
29	1400	1050	156.000	82.278	144	104	74.867	4	3	N	P	CVT	○	○	○
30	1600	1200	161.000	74.537	168	112	59.869	4	3	N	P	CVT	○	○	○
31	1600	1200	162.000	75.000	192	64	60.000	3	1	P	P	VESA Standard	○	○	○
32	1600	1200	175.500	81.250	192	64	65.000	3	1	P	P	VESA Standard	○	○	×
33	1600	1200	189.000	87.500	192	64	70.000	3	1	P	P	VESA Standard	○	○	×
34	1920	1080	141.500	55.621	200	112	49.929	5	3	N	P	CVT	○	○	○
35	1920	1080	138.500	66.587	32	48	59.934	5	2	P	N	CVT Red.Blanking	○	○	○
36	1920	1200	154.000	74.038	32	48	59.950	6	3	P	N	CVT Red.Blanking	○	○	○

●VIDEO Signal (CVBS、S-VIDEO)


No.	Resolution		Dot CLK (MHz)	Horizontal frequency			Vertical frequency			Polarity		Signal name Type	CVBS S-VIDEO
	Horizontal	Vertical		H Freq	Sync Time	Front Porch	V Freq	Sync Time	Front Porch	H	V		
	(Dot)	(H)		(KHz)	(Dot)	(Dot)	(Hz)	(H)	(H)				
1				15.734			59.940					NTSC	○
2				15.625			50.000					PAL	○

When signals other than the listed above are input, the monitor does not will not display properly or noise appears. It is not a malfunction of the monitor. Please input the standard signals listed above.

Caution : When inputting unique signals other than listed above, please contact our company's sales department.

* The information in this documents is subject to change without notice.

VI. External Dimensions

* The information in this documents is subject to change without notice.

VII. Attachments

1. Nylon clamps (NK-3)	1 pieces
2. Bind screws (Φ38mm)	1 pieces
3. AC Adaptor	1 unit
4. Instruction Manual	1 set
5. Quick Reference Guide	1 set

VIII. Warranty

- The warranty period for this product is 1 year.
- After repair service, warranty period for replaced parts are 6 months.
- As a rule, I will assume seven years after the final production of products is completed.
However, when a key part(LCD panel,etc.) was stopped producing, it is inapplicable.
- In principle, malfunction products shall be send back to our company's factory for repair.
The customers shall bear the cost (business trip expense and technology fee of the monitor detaching, etc.). Sender pays for the product costs.
- This warranty does not apply ;
 - : To damages caused by misuse and authorized repairs and modifications.
 - : Malfunction and damaged caused by using other companies' products other than our company's products.
 - : Malfunction and damaged caused by using other maintenance parts other than our company designated original factory parts to carry out maintenance and repair.
 - : Malfunction and damaged caused by not following caution items and operation methods described in the user Instruction Manual and Quick Reference Guide that is attached to the product.
 - : Malfunction and damaged caused by using power supply and installation environment etc. that do not meet this product required use conditions as described in user Instruction Manual and Quick Reference Guide.
 - : To damages caused by accident, abuse, flood, fire, earthquake, other external caves, salt age, gas harm, or abnormal voltage.
 - : To damages caused by the product being dropped, transported, etc.

After the warranty period, the warranty period for the repair of functional parts etc. questions, please contact the store you purchase this product or our company's service department. The parts that may be continuously used after repair will be charged and repaired.

●Warranty parts

Please contact the store you purchased your products or our company's service department for the basis of replacement warranty parts.

1) Backlight

Caution: Backlight replacement work shall be carried out by panel maker.

2) Exclusive use AC Adaptor

Caution

- No part of this document may be reproduced in any means for any purpose.
 - The information in this documents is subject to change without notice.
-

Cases For Indemnity (Limited Warranty)

We shall be exempted from taking responsibility and held harmless for damage or losses incurred by the user in the following cases.

- In the case damage or losses are caused by fire, earthquake, or other acts of God, acts by a third party, deliberate or accidental misuse by the user, or use under extreme operating conditions.
 - In the case of indirect, additional, consequential damages (loss of business interests, suspension of business activities) are incurred as result of malfunction or non-function of the equipment, we shall be exempted from responsibility for such damages.
 - In the case infringement of third party intellectual properties by use of this unit, we have been exempt from the responsibility for such infringement except directly related to our company's product configuration and production method.
 - In the case damage or losses are caused by failure to observe the information contained in the instructions in this Instruction Manual and Quick Reference Guide and specifications.
 - In the case damage or losses are caused by use contrary to the instructions in this Instruction Manual and Quick Reference Guide and specifications.
 - In the case damage or losses are caused by malfunction due to connection to other equipment.
 - In the case damage or losses are caused by malfunction or other problems resulting from use of equipment or software that is not specified.
 - In the case damage or losses are caused by repair or modification conducted by the customer or any unauthorized third party (such as an unauthorized service representative).
 - Expenses we bear on this product shall be limited to the individual price of the product.
 - The warranty is limited to the item described in this use instruction.
-

Restrictions for use

- Should the equipment be used in the following conditions or environments, give consideration to safety measures and inform us of such usage:
 1. Use of the equipment in the conditions or environment contrary to those specified, or use outdoors.
 2. Use of the equipment in applications expected to cause potential hazard to people or property, which require special safety measures to be adopted.
- Determination of applicability of equipment or devices concerned shall be determined after analysis or testing as necessary by the designer of such equipment or devices, or personnel related to the specifications. Such designer or personnel shall assure the performance and safety of the equipment or devices.
- This product is not designed or manufactured to be used for control of equipment directly concerned with human life (※1) or equipment relating to maintenance of public services/functions involving factors of safety (※2). Therefore, the product shall not be used for such applications.

Accordingly, we disclaims any all liability arising out of the use of the products in the listed above applications.

(※1): Equipment directly concerned with human life refers to.

- Exhaust control equipment for exhaust gases such as toxic fumes or smoke.
- Equipment mandatory to be installed by various laws and regulations such as the Fire Act or Building Standard Law
- Equipment related to the above

(※2): Equipment relating to maintenance of public services/functions involving factors of safety refers to.

- Equipment for nuclear power generation
- Equipment related to the above

* The information in this documents is subject to change without notice.

Request in the use

- The handling carefully

Not fall down or collide. It may lead to machine broken down and malfunction.

Excessively do not bend or apply force to the cable. It may broke the cable or cause malfunction of the unit.

- Ambient temperature and ambient humidity

Please do not use it at the place of ambient temperature and humidity more than the specifications of this product. Picture quality deteriorates and it is bad influence in an internal part. In addition, please be careful to the use at a place places with direct sunlight.

- Picture quality

The image may show unevenness. This is due to the characteristic of LCD and not malfunction.

The image may be flickering when displaying checkered pattern with every dot in half tone or at middle luminance level. This is not malfunction of the monitor.

The image may show red, green or dots, or red, green or blue non-flashing dots. This is not malfunction of the monitor.

- Burn out

When the same picture is displayed for a long period of time, the lag of the previous picture may appear during changing images. The lag will be disappeared gradually. However the screen may be burn out and the burn out image stay on the screen permanently when the same picture has been displayed for extremely long period. Please be careful for screen the burn out.

Publication

Version 2
Specifications

completed 2008-11
T24MSA001-MD

TOSHIBA TELI CORPORATION

Head Office : 7-1, 4 chome, Asahigaoka, Hino-shi, Tokyo, 191-0065, Japan
(Overseas Sales Department)

Phone : +81-42-589-8771

Fax : +81-42-589-8774

URL : <http://www.toshiba-teli.co.jp>